

ORDER FORM



Invoicing Details	Shipping Details					
Company name:	Ship to:					
Delivery date required:/						
Delivery Address:	CUSTOMER STAMP AND SIGNATURE					
P.Code: Town:						
Tel.: Fax:						
# Order: Date of the Order:/						

Pro Medicare srl grants 24 months warranty on production defects from the date of the first put into service and 12 months on covers and components subjected to wearing. Warranty is void if the product is used in any other way than stated in the Instruction Manual.

Note: Any request for variation in dimensions, materials and type compared to the standard does not have CE marking, therefore the professional user becomes the manufacturer and has the obligation to draw up the technical documentation of the device. Images are purely illustrative and may not fully reflect reality.

Standard Configuration

- Frame Adacta Liga Light anodised satin
- · Fixed folding down back posts
- Reclining mechanism by predetermined plate w/pin 3°,9°,15°,22°
- Continuous Tilt in Space by drive lever from 2° up to 35° (from 6° up to 39°)1 $1 \\ \textit{Tilt in space can be up-graded of 4}^{\circ} \textit{ with an increment of the starting position of the frame by operating on the positioning}$ holes of the caster fork
- Adjustment of the seat depth
- Height adjustment: High (49 cm from the floor)/low (43 cm from the floor)* *(with 300mm and 500mm wheels only)
- Detachable footplates with 90° fixed knee angle
- Fixed Anti-tip system
- Wheels Kit with brake
- Rigid wooden seat base covered leather flex
- Multi-Adjustable padded hip quides
- Multi-Adjustable padded arm supports
- Versa Inserto Trunk/Pelvis positioning system
- Capitis Headrest
- Fixatis Pelvic Belt



reflect the standard configuration.

□ = Optional with up-charge **◊= Standard** O = Optional without up-charge

* If a selection is not made, the product will be delivered in its standard configuration Frame Color Frame Size (Outside tubes Width x Depth cm) Frame Type

O Turquois O Pink	V RIGIA	O US (30x38)	O XXS (35x38)	O XS (35x42)	O XS0 (38x42)					
V Tarquois OTHIK	O Folding	0 03 (30x30)	0 XX3 (33X30)	O N3 (33X42)	0 X30 (30X+2)					
□ Pre-assembly of postural components (collars for flat base, backrest attachments, headrest hardware)										
Knee Angle	◊ 90° ○ 75° □ elevating									
Castors	◊ 150mm (6") poly grey	O black								
Castors	O 175mm (7") poly O grey	O black								
Posterior Wheels (height 43 cm from seat to the floor with 300 mm and 500 mm wheels; height 47 cm from seat to the floor with 600 mm wheels)										

§ 300mm (12") poly grey with parking brake lever

- O 300mm (12") poly black with parking brake lever
- O 300mm (12") pneumatic grey with parking brake lever
- □ 300mm (12") poly black with drum brake
- □ 300mm (12") poly grey with drum brake
- □ 300mm (12") pneumatic grey with drum brake
- □ 300mm (12") poly black with drum brake and parking brake lever
- □ 300mm (12") poly grey with drum brake and parking brake lever
- □ 300mm (12") pneumatic grey with drum brake and parking brake lever
- □ 400mm (16") poly black with parking brake lever (only 45 cm height)
- □ 400mm (16") poly black with drum brake (only 45 cm height)
- □ 400mm (16") poly black with drum brake and parking brake lever (only 45 cm height)
- □ 500mm (20") poly grey without handrim with parking brake lever
- □ 500mm (20") pneumatic grey without handrim with parking brake lever

□ handrim for 500mm (20") wheels

- □ 500mm (20") pneumatic grey with handrim with drum brake
- □ 500mm (20") pneumatic grey with handrim with drum brake and parking brake lever
- □ 600mm (24") pneumatic grey with handrim with parking brake lever (not available on size US)

(Note: The tilt in space is reduced to 20°, height from the floor 47 cm)

Mod.Ord.Ver.Tit. S.P EN Rev.18 10/2025 1/3

□ Inserto Seat Domino										
□ Inserto Seat Domino										
□ Inserto Seat Domino										
□ Inserto Seat Domino										
□ Inserto Seat Domino										
Back Supporto										
 □ Kit for telescopic and angle adjustable single handlebar □ Kit for angle adjustable handlebar for folding frame 										

Accessories

Thoracic Supports for Inserto Back Supporto *Please, indicate the type, the shapes and the size desired

moracie supports for miscres sack sup	Trease, maleate the type, the shapes and the	ic size desired	_					
Type of Bracket	Shape of the Pad			cm	☐ size 1	☐ size 2	☐ size 3	☐ size 4
Fixed	Flat			н	6	8	10	12
□ RH □ LH	□ RH □ LH			L	6,5	8,5	10,5	13,5
Swing-Away	Contoured			cm	☐ size 1	☐ size 2	☐ size 3	
□ RH □ LH	□ RH □ LH			н	6	8	10	
			-	L	10,5	12,5	14,5	

Thoracic Supports for Inserto Back Effecto Lato *Please, indicate the type, the shape of the brackets, the pad and the size desired

The fact Supports for inserted back Effects Eath Fredse, indicate the type, the shape of the brackets, the pad and the size desired										
Type of Bracket		Shape of the Pad		cm	☐ size 1	☐ size 2	☐ size 3	☐ size 4		
Fixed		Flat			н	6	8	10	12	
□ RH □ LH		□ RH □ LH			L	8	10	12	15	
Swing-Away		Contoured	·MD		cm	☐ size 1	☐ size 2	☐ size 3		
□ RH □ LH		□ RH □ LH			н	6	8	10		
Detachable		Bracket shape (for the selection see TECHNI	CAL INFORMATION)		L	12	14	16		
□ dx □ sx		type A: □ RH □ LH type B:	□ RH □ LH			-	<u> </u>		•'	

Pedal for barriers	□ RH □ LH		Abductor w	redge on	ly with rig	id seat bo	ise		
☐ Vehicle Transport		cm	□S	□М	□ removable with				
' ''	wiith 300mm and 500mm poster	ior wheels) heelchair as a seat in the vehicle if the		Н	7.5	9	lever		
weight of the user is lower th		leerchair as a seat in the vehicle ij the		D	9	10	□ removable with		
Lower limb accesso	ries			L	7.5	9	mounting block		
☐ Separate Calf Sup	ports swing-away		Note: if not ind	icate the a	bductor we	dge will be	removable with lever		
☐ Padded Calf Strap	•		Please indic	ate the i	position	of holes	under the seat		
Cover for legrest tu	be □ RH □ LH		(if not indicated the abductor wedge will be fixed in the standard position)						
Cover for Knee Ang	gle joint □ RH □ LH	ł	□ standard						
☐ Padding and cove	er of footplate in leath	herette	not centred □ on the RH □ on the LH						
Capitis Headrest Co	mfort size:		Tray						
O Small	◊ Medium	O Large	□ Rectangu	lar Woo	den Tray	/ w/late	ral flip hardware		
Fixatis Pelvic belt 2	☐ Rectangular Wooden Tray w/detachable hardware								
o xs	◊ <i>S</i>	0 M	□ Padding	for Tray					
Fixatis Dynamic But	terfly Harness				•	•			
□XS	⊓S	⊓М	1						

Mod.Ord.Ver.Tit. S.P EN Rev.18 10/2025 2/3

☐ Kit of ABS laterals multi-adjustable, detachable and w/multi-adjustable arm support

in width and depth with 2.5 cm gap

make a selection:

o small arm support (from US to XSO) (dimensions: 5 cm x 25 cm)

o large arm support (from XS1 to XL) (dimensions: 6 cm x 30 cm)

make a selection:

o short vertical bar (from 22 cm to 31 cm*)

o long vertical bar (from 29 cm to 38 cm*)

* distance from armrest to the seat rail



☐ Multi-adjustable and detachable hips and arms positioning kit

- pelvis (width, depth, height, adduction/abduction adjustable)
- arm support (width, depth, height, inclination, intra/extra-rotation adjustable)

For independent adjustments referring to the width the gap is 3 cm

The armrests and the hip guides are assembled on the frames as the following indications without possibilities to be modified

(armrests dimensions: 8 cm x 22 cm) o **small** (from frame size US up to XXS) (hips guide dimensions: 25 cm x 10 cm with slot of 5 cm)

○ <u>large</u> (from frame size XS to XSO) (hips guide dimensions: 30 cm x 12 cm with slot of 7.5 cm) (armrests dimensions: 10 cm x 30 cm)



Sizes of the arm supports and hip guides. Note: The components cannot be interchangeable

Sizes of the arm supports a	illu liip g	ulues. IV	ote. The	compon	ents cum					
Linear Hip Guides										
	cm	US	XXS	XS	XS0					
"	Н	10	10	10	10					
s= 2,5 cm	L	25	25	30	30					
Flat Arm Support										
s=2,5	cm	US	XXS	XS	XS0					
	Н	10	10	10	10					
4	L	25	25	30	30					

gea.are											
ABS Lateral	•										
1	cm	US	XXS	XS	XS0						
H	Н	11	11	11	13						
L L	L	24	24	24	28						
Armrest of the ABS Late	Armrest of the ABS Lateral kit										
s=2,5	cm	US	XXS	XS	XS0						
	Н	5	5	5	5						
H	L	27	27	27	33						

	Technical Information													
Ti	itti Size	Maximun Load (kg)	Frame Weight (kg)	Pivot Width	Maximum overall dimensions in width/length (90° knee angle and footrest positioned forwards) Values expressed in cm									
		(1/6)	(1.6/	(cm)	n) standard 300mm 300mm wheels standard 400mm 400mm wheels standard 500mm 500mm wheels						600mm wheels			
US	Rigid	50	14.1	69	48/96	49/96	50/99	53/99	50/101	53/101	,			
03	Folding	50	14.9	09	46/90	49/90	30/99	25/33	30/101	55/101	/			
xxs	Rigid	50	14.4	69	53/96	54/96	55/99	58/99	55/101	58/101	55/104			
^^3	Folding	50	15.3	09	33/30	34/90	55/55	36/33	55/101	36/101	55/104			
xs	Rigid	50	14.7	72	53/99	54/99	55/100	58/100	55/102	58/102	55/107			
^3	Folding	50	15.6	/2	33/99	34/99	33/100	36/100	55/102	58/102	55/10/			
XSO	Rigid	50	15.4	75	56/99	57/99	58/100	61/100	58/102	61/102	58/107			
^30	Folding	50	15.9	/5	30/99	37/99	36/100	01/100	36/102	01/102	36/10/			

Note: Frame weight only with wheels, knee angle and footrest as standard configuration. The length dimensions, referred to the average values, are subjected to variations due to the displacement of the wheel plates. The values reported are intended for the standard configuration.

	User Pai	rameters	Frame Depth	Cushion Size	Supporto Back sizes (cm)			Effecto Lato Back sizes (cm)				
Titti Size	Pelvis Width (cm)	Seat Depth (cm)	Adjustments (cm)	(IW x D) (cm)	Backrest	Height	Shell Width	Trunk Width	Backrest	Height	Shell Width	Trunk Width
US	20-30	28-38	da 36 a 44	30x38	3033	33	27	17-23	3033	33	30	w/bracket A 18-28 w/bracket B 13-23
xxs	26-36	28-38	da 36 a 44	36x38	3336	36	30	20-26	3338	38	33	w/bracket A 21-31 w/bracket B 16-26
xs	26-36	32-42	da 40 a 48	36x42	3640	40	33	23-29	3642	42	36	w/bracket A 24-34 w/bracket B 19-29
XSO	29-39	34-44	da 40 a 48	39X45	3944	44	36	26-34	3945	45	39	w/bracket A 27-37 w/bracket B 22-32

Pro Medicare)

Pro Medicare s.r.l.

Via Montagna Z.I. Lotto 41 - 72023 Mesagne (Br) - ITALY Website: www.promedicare.it - E-mail: info@promedicare.it: sales@promedicare.it

Customer Service **(+ 39) 0831. 777 840**

Mod.Ord.Ver.Tit. S.P EN Rev.18 10/2025 3/3